

SI - 1000

SI - 1000

SI - 1000

The SI-1000 has the same high quality as the Legacy-500 but is designed for the clinic where time is a factor. With an internal timer you never worry time will vary, a simple push button starts table to desired run time every time.



Standard Features:

- 24' wide x 74' long
- 10-minute internal timer (set at desired run time)
- Push button start
- Counter-clock massage direction
- Adjustable massage pressure
- Automatic massage pressure lowering
- Treat full spine
- Cervical roll and knee bolster

Options:

- Reverse side controls
- Therapy top
- 27' over sized top

Specifications:

- Electrical: 120V, 60Hz, 10 Amp.
- Weight Capacity: 400 lbs.
- One year limited warranty
- Life time warranty on rollers
- Shipping weight 349 lbs.

SI-1000 Table	\$2,375.00
RSC Reverse Side Control	\$250.00
ADM Admiral Upholstrey	\$300.00
27" Over Sized Top	\$150.00
Patient Straps and Rails	\$250.00
Spring Loaded Roller Pressure	\$350.00

Legacy-500 Intersegmental Traction Table

Legacy-500 Intersegmental Traction Table

Legacy-500 Intersegmental Traction Table

With 2" more travel and almost 2" wider rollers than any one in the industry the Legacy-500 is specially designed to release muscle tension and deliver a soothing massage. Under the thick, soft, padded top, three large rotating rollers travel the back to release tension and relax tired muscles.



Standard Features:

- 24" wide x 74" long
- 30-minute timer
- Reversible massage direction
- Adjustable massage pressure
- Vibration with On-Off feature
- Automatic massage pressure lowering
- Treat specific area or full spine
- Cervical roll and knee bolster

Options:

- Reverse side controls
- Therapy top
- 27" over sized top
- Patient Straps and Rails
- Spring Loaded Roller Pressure
- Push Button Timer, 1-10 Minute Range
- Patient Hand Control - Stops Travel and control Roller Pressure

Specifications:

- Electrical: 120V, 60Hz, 10 Amp.
- Weight Capacity: 400 lbs.
- One year limited warranty
- Life time warranty on rollers
- Shipping Weight 349 lbs.

Legacy-500 Intersegmental Traction Table 25" high, 24" wide, 74" long, includes 8" wide massage rollers, 30-minute timer, Vibration, adjustable massage pressure, Reversible roller travel, isolated or full spine massage, automatic massage lowering device, cervical gel roll and knee bolster, 400 lb. capacity, Lifetime warranty on rollers, One year limited warranty.	\$2,500.00
Legacy-501 Intersegmental Traction Table 25" high, 24" wide, 74" long, includes 8" wide massage rollers, 30-minute timer, Vibration, adjustable massage pressure, Reversible roller travel, isolated or full spine massage, automatic massage lowering device, cervical gel roll and knee bolster, 12" face opening with face cushions, 400 lb. capacity, Lifetime warranty on rollers, One year limited warranty.	\$2,600.00
Legacy-500 and Legacy-501 Intersegmental Traction Table Options	
LEG-RSC Reverse side control	\$250.00
LEG-ADM Admiral upholstery	\$300.00
27" over sized top	\$150.00
Patient Straps and Rails	\$250.00
Spring Loaded Roller Pressure	\$350.00
Push Button Timer, 1-10 minute range	\$300.00
Patient Hand Control	\$300.00

Color Chart

/ Color Chart

Upholstery Options

Color is not exact.

	Vinyl Upholstery		Admiral Upholstery	
837 Oxblood		833 Crimson		A57 Colonial Blue
844 Flint Gray		1009 Mauve		A83 Burgundy
830 Walnut		832 Blue		A30 Chocolate
9375 Charcoal		1011 Jadestone		A10 Charcoal
1018 Blue Heron		835 Black		A12 Platinum Gray
839 Pecan				
8317 Mocha Brown				

MUST specify at time of order, NO changes allowed