

Ordering Information

Item#	Description
Intellect TranSport® Ultrasound	
2782	Intellect Transport Ultrasound
Standard Accessories	
27383	5 cm ² Sound Head Applicator
4248	Conductor™ Ultrasound Gel 8.5 oz (250 ml) Bottle (1x)
Optional Accessories	
27478	Intellect Transport Nickel Metal Hydride Battery Pack
27381	1 cm ² Sound Head Applicator
27382	2 cm ² Sound Head Applicator
27384	10 cm ² Sound Head Applicator
27467	Intellect Transport Carry Bag
2780ASY	Therapy Cart (assembled)
2884	Cart Adapter (Gray) Allows the Intellect Transport unit to attach securely to the Therapy Cart.

Note: The Cart Adapter, part number 2884, is necessary when ordering a cart for a Transport Unit.

Technical Specifications

Intellect Transport Ultrasound

Mains Power: 100-240 VAC-50/60 Hz, 1.0 A

Weight: 5 lb (2.3 kg) /Shipping - 9 lb (4 kg)

Dimensions*: 12.9" x 11.3" x 6.4" (32.8 x 28.8 x 16.3 cm)

Electrical Safety Class: Class 1, Type B

Safety Tests: UL/IEC/EN 60601-1, IEC/EN 60601-1-2, IEC 60601-2-5

* Dimensions are listed as L x W x H

Intellect TranSport® Ultrasound



Features

Intellect TranSport® Ultrasound

- Fully functional 1 and 3.3 MHz frequencies
- Ten user-defined memory positions for User Protocols
- Pulsed and continuous therapy operation (10%, 20%, 50% and 100%)
- Ergonomically designed Sound Head Applicators in variety of sizes (1 cm², 2 cm², 5 cm² and 10 cm²)
- Sound heads are watertight for use in underwater therapy
- All Sound Head Applicators are interchangeable due to Electronic Signature™ residing on all transducers
- Head-warming feature on all Sound Head Applicators
- Unit designed for table mount, therapy cart or mobile use
- 110V or battery powered option
- Optional Therapy Cart integrates with the unit to help maximize storage of clinical supplies
- Two-year warranty



Cart Adapter



Optional Battery Pack



Intellect TranSport® Electrotherapy



Features

Intellect TranSport® Electrotherapy

- Two channels of electrical stimulation output
- Logical control system and a large easy-to-read graphical LCD display
- Four standard waveforms on board: Interferential, Premodulated, High Volt, and Russian
- Ten user-defined memory positions for User Protocols
- Independent intensity and parameter controls for each channel
- Unit designed for tabletop, wall mount, therapy cart or mobile use
- 110V or battery-powered option
- Optional Therapy Cart integrates with the unit to maximize storage of clinical supplies
- Two-year warranty



The custom designed carrying bag (optional) is lightweight. Protect and transport the unit and accessories for therapy on the go.

Ordering Information

Item#	Description
Intellect TranSport® Electrotherapy	
2783	Intellect Transport® 2-Channel Electrotherapy
Standard Accessories	
27312	Electrotherapy Leadwire Channel 1
27313	Electrotherapy Leadwire Channel 2
Optional Accessories	
27478	Intellect Transport Nickel Metal Hydride Battery Pack
72852	Black Rubber Carbon 3" (8 cm) Electrodes (2x)
72853	Red Rubber Carbon 3" (8 cm) Electrodes (2x)
10648	Nylatex® Wrap 2.5" x 24" (8x)
27467	Intellect Transport Carry Bag
2780ASY	Therapy Cart (assembled)
2884	Cart Adapter (Gray) Allows the Intellect Transport unit to attach securely to the Therapy Cart.
Note: The Cart Adapter, part number 2884, is necessary when ordering a cart for a Transport Unit.	

Waveforms

	INTERFERENTIAL <ul style="list-style-type: none"> • Continuous • Manual Vector Scan
	PREMODULATED <ul style="list-style-type: none"> • Continuous • Cycle Timed
	HIGH VOLT <ul style="list-style-type: none"> • Twin peak
	RUSSIAN <ul style="list-style-type: none"> • Continuous • Cycle Timed

Technical Specifications

Intellect Transport Electrotherapy	
Mains Power:	100-240 VAC-50/60 Hz, 1.0 A
Weight:	5 lb (2.3 kg) /Shipping - 9 lb (4 kg)
Dimensions*:	12.9" x 11.3" x 6.4" (32.8 x 28.8 x 16.3 cm)
Electrical Safety Class:	Class 1, Type B
Safety Tests:	UL/IEC/EN 60601-1, IEC/EN 60601-1-2, IEC 60601-2-5

* Dimensions are listed as L x W x H

Ordering Information

Item# Description

Intelect TranSport® Combo

2738 Intelect Transport Combo

Standard Accessories

27383 5 cm² Sound Head Applicator

4248 Conductor™ Ultrasound Gel 8.5 oz (250 ml) Bottle (1x)

27312 Electrotherapy Leadwire Channel 1

27313 Electrotherapy Leadwire Channel 2

Optional Accessories

27478 Intelect Transport Nickel Metal Hydride Battery Pack

72852 Black Rubber Carbon 3" (8 cm) Electrodes (2x)

72853 Red Rubber Carbon 3" (8 cm) Electrodes (2x)

10648 Nylatex® Wrap 2.5" x 24" (2x)

27381 1 cm² Sound Head Applicator

27381 2 cm² Sound Head Applicator

27384 10 cm² Sound Head Applicator

27467 Intelect Transport Carry Bag

Waveforms



INTERFERENTIAL

- Continuous
- Manual Vector Scan



PREMODULATED

- Continuous
- Cycle Timed



HIGH VOLT

- Twin peak



RUSSIAN

- Continuous
- Cycle Timed

Technical Specifications

Intelect TranSport® Combo

Mains Power: 100-240 VAC-50/60 Hz, 1.0 A

Weight: 5 lb (2.3 kg) / Shipping - 9 lb (4 kg)

Dimensions*: 12.9" x 11.3" x 6.4" (32.8 x 28.8 x 16.3 cm)

Electrical Safety Class: Class 1, Type B

Safety Tests: UL/IEC/EN 60601-1, IEC/EN 60601-1-2, IEC 60601-2-5

* Dimensions are listed as L x W x H

Intelect TranSport® Combo



Features

Intelect TranSport® Combo

- Fully functional 1 and 3.3 MHz ultrasound frequencies
- Two channels of electrical stimulation output
- Ten user-defined memory positions for User Protocols
- Four standard waveforms on board: Interferential, Premodulated, High Volt, and Russian
- Pulsed and continuous therapy operation (10%, 20%, 50% and 100%)
- Independent intensity and parameter controls for each channel
- Ergonomically designed Sound Head Applicators in a variety of sizes (1 cm², 2 cm², 5 cm² and 10 cm²)
- All Sound Head Applicators are interchangeable with Electronic Signature™ residing on all transducers
- Watertight sound heads for use in underwater therapy
- Head-warming feature on all Sound Head Applicators
- Unit designed for tabletop, wall mount, therapy cart or mobile use
- 110V or battery powered option
- Two-year warranty

