



402 Flexion Stationary



Table Specifications

Height Range 18 inches thru 30 inches

Length Extended 78in.

Length Retracted 67in.

Cushion Width 21in.

Table Width 22in.

Shipping Weight
Approx. 200lbs.

Actual Weight
Approx. 150lbs.

The 402 Flexion-Stationary

table features a narrowed head and shoulder section for added patient comfort. A simple control lever provides instant locking of flexion or lateral bending. Hydraulic control provides more constant resistance through the flexion range, regardless of the patient's size and weight. **Lloyd's drop** design allows the amount of force applied to be equal, whether that force is applied on the right, left or center of the pelvic cushion, making them superior to other drops on the market. The automatic cocking drops are built into the table with controls on each side and can be set manually in the event of electric failure. No air compressor is required. A specially designed spring release on drop chest and lumbar and flexion T-BAR that retracts downward or sideways are optional.

- Tilting headpiece with add-on paper attachment
- One piece adjustable chest
- Manual Flexion Distraction
- Adjustable Ankle Rest (Horizontal)
- Axial extension & retraction of pelvic section (Manual Foot Pump)
- Lateral bending
- Ankle Straps
- **Optional features**
- Lloyd's superior designed manual cocking drops
- Manual drop Ankle
- Pelvic rotation
- Total Dimension Flexion Headpiece
- Grab Bar
- Spring Breakaway Chest
- Rotating chest lumbar with breakaway
- Retractable flexion T-BAR
- Angled T-Bar
- Snap-on Covers

