



Galaxy Elevation

Table Specifications

Height Range	19" to 36"
Fully Extended	83" long
Fully Retracted	62" long
Cushion Width	21"
Table Width	22"
Shipping Weight	400 lbs. (approx)
Actual Weight	350 lbs. (approx)



The Galaxy Elevation raises the patient from 19" to 36" virtually eliminating stooping and bending. Doctors using this table find they are less fatigued at the end of the day. Optional automatic and manual drops are available on all sections even ankle! Fully adjustable headrest and elevating, extendable ankle rest position manually to accommodate patients 5' to 7'.

- Hydraulically Controlled Elevation
- Deluxe Adjustable Headpiece
- One piece adjustable chest
- Spring Breakaway Chest
- Adjustable Chest (vertical)
- Adjustable Pelvic (vertical)
- Adjustable Ankle (vertical & horizontal)

Optional Features

- Lloyd's superior designed manual & auto cocking drops
- Manual & auto drop Ankle
- Total Dimension Flexion Headpiece
- Electronic Height Select
- Snap-on Covers



Imaging Solutions

1-844-944-6243

MIKE@TRU-DIGITAL-IMAGING.COM