



Galaxy Hylo

Hydraulically Controlled Hylo Table




Table Specifications

Height Range	18" to 28"
Fully Extended	84" long
Fully Retracted	76" long
Cushion Width	21" long
Table Width	22" long
Shipping Weight	400 lbs. (approx)
Actual Weight	350 lbs. (approx)

The Favorite Workhorse of the Chiropractic Profession. Standard features include hydraulic tilt from horizontal to 70 degrees, infinitely adjustable speed and hydraulic extension retraction. A chest lumbar drop with spring breakaway, Pelvic Drop, ankle drop, plus many more options are available. Also available on the Galaxy Hylo is the All-New Lloyd Total Response Drop system on all sections. The Lloyd Total Response Drop System features a one-piece chest cushion while offering a dorsal and lumbar drop that cock and drop individually.

- Hydraulically Controlled Hylo
- Deluxe Adjustable Headpiece
- One piece adjustable chest
- Extension / Retraction of Length (Power)
- Pelvic Separation
- Adjustable Chest (vertical)
- Adjustable Pelvic (vertical)
- Tilt Safety Switch
- Speed Adjustment for Hylo feature
- Breakaway Chest

Optional Features

- Lloyd's superior designed auto and manual cocking drops
- Manual & Automatic drop Ankle
- Total Dimension Flexion Headpiece
- Rotating chest lumbar with breakaway
- Multi-Directional Drop Pelvic
- Snap-on Covers



TRU-DIGITAL

Imaging Solutions

1-844-944-6243

MIKE@TRU-DIGITAL-IMAGING.COM