



# Galaxy McManis Elevation

## Hydraulically Controlled Elevation Table



### Table Specifications

Fully Extended 90 x 22in

Fully Retracted 64 x 22in

Fully Raised Ht. 24in. OR 40in.

Fully Lowered Ht. 22in. OR 38in.

Ship Wt. Approx 430lbs

Actual Wt. Approx. 380lbs

Unequaled versatility in an all-purpose flexion table. Manual flexion and powered pelvic extension and retraction are standard. Optional features include: electronic height selector, and programmable automatic flexion, rotating chest, rotating pelvic, manual & automatic cocking drops on all sections including the ankle!

**Lloyd's drop** design allows the amount of force applied to be equal, whether that force is applied on the right, left or center of the pelvic cushion, making them superior to other drops on the market. The automatic cocking drops are built into the table with controls on each side and can be set manually in the event of electric failure. No air compressor is required. A specially designed spring release on drop chest and lumbar and flexion T-BAR that retracts downward or sideways are optional. Optional manual and automatic drops are available on an optional rotating chest lumbar section.

- Hydraulically controlled Elevation
- Deluxe / Adjustable Headpiece
- One piece adjustable chest
- Chest & Pelvic Separation
- Extension / Retraction of Length
- Semi-automatic and manual flexion distraction
- Adjustable Ankle Rest (Vertical and Horizontal)
- Powered axial distraction
- Lateral bending
- Ankle Straps

### Optional features

- Automatic Programmable Flexion
- Lloyd's superior designed auto and manual cocking drops
- Electronic Height Select
- Pelvic rotation
- Total Dimension Flexion Headpiece
- Grab Bar
- Rotating chest lumbar with breakaway
- Retractable flexion T-BAR
- Angled T-BAR
- Span-on covers

