



# Galaxy McManis HYLO

*A Safety Oriented Flexion Table That Stresses Comfort*

## Table Specifications

Height Range 22" to 28"

Fully Extended 91in.

Fully Retracted 65in.  
including footplate

Fully Retracted 65in.  
excluding footplate

Cushion Width 21in.

Table Width 22in.

Ship Wt. Approx 400lbs

Actual Wt. Approx. 350lbs



The Galaxy McManis Hylo is a flexion distraction table designed with both doctor and patient in mind. Available in 22" to 28" heights, it features fully adjustable head, chest, lumbar and ankle sections. Simple control levers provide precision control of the pelvic flexion with lateral bending and traction.

The standard chest lumbar section, with spring breakaway, is designed to comfortably accommodate pregnant or large abdominal patients. A lumbar extension combined with extendable and retractable pelvic and ankle sections makes it possible to accommodate patients from 5' to 6' 11".

- Hydraulically controlled Hylo (tilt)
- Deluxe / Adjustable Headpiece
- One piece adjustable chest
- Chest & Pelvic Separation (Power)
- Extension / Retraction of Length (Power)
- Semi-automatic and manual flexion distraction
- Adjustable Ankle Rest (Vertical and Horizontal)
- Tilt Safety Switch
- Speed adjustment of Hylo Feature
- Powered axial distraction
- Lateral bending
- Ankle Straps

## Optional Features

- Lloyd's superior designed auto and manual cocking drops
- Manual & Automatic drop Ankle
- Pelvic rotation
- Automatic Flexion Distraction
- Total Dimension Flexion Headpiece
- Grab Bar
- Rotating chest lumbar with breakaway
- Retractable flexion T-BAR
- Angled T-BAR
- Snap-on covers
- Manual & Automatic drop ankle

