



Galaxy Ultimate Table

Elevation | HYLO | Flexion Distraction

Table Specifications

Height Range 22" to 38"

Fully Extended 91in.

Fully Retracted 81in.
including footplate

Fully Retracted 65in.
excluding footplate

Cushion Width 21in.

Table Width 22in.

Shipping Wt. Approx
500lbs

Actual Wt. Approx. 450lbs



The Galaxy Ultimate can be smoothly tilted or lowered hydraulically at infinitely adjustable speeds from horizontal to 70° controlled by foot pedals on both sides of the table. A safety switch prevents the table from accidentally lowering from the hylo position. A foot plate automatically comes up to meet the patient's feet when table begins to lift.

- Hydraulically controlled elevation
- Tilting (Hylo) capability
- Automatic, semi-automatic manual flexion distraction
- Electronic height selector
- Powered axial distraction
- Lateral bending
- One piece adjustable chest with spring breakaway
- Extension/Retraction of length - powered
- Electric height selector
- Deluxe/Adjustable Headpiece
- Chest/Pelvic separation (powered)
- Adjustable ankle rest - vertical & horizontal
- Pelvic - Flexion/Distraction (manual, semi-automatic & automatic)
- Pelvic - Lateral side bending
- Axial Extension/Retraction (powered)
- Ankle straps
- Tilt safety switch
- Speed adjustment for hylo feature
- Optional Features
 - Lloyd's superior designed auto and manual cocking drops
 - Pelvic rotation
 - Rotating chest lumbar with breakaway
 - Retractable flexion T-BAR
 - Snap-on covers
 - Total dimension flexion headpiece
 - Grab bar
 - Angled T-bar



1-844-944-6243
MIKE@TRU-DIGITAL-IMAGING.COM