



DC-2001 Table



Table Specifications

Electrical: 120V, 60Hz, 10amp

Weight Capacity 300lbs

Warranty 1 year limited

Shipping Weight 250lbs

Table Dimensions

Length 65in. to 73in. with tilting head

Length 59in to 68in with elevating head

The DC-2001 is designed with both you and the patient in mind. The DC-2001 adjusts from 21" to 29" or the optional height range of 18" to 28". You can customize your DC-2001 from the many available options.

Standard Features

- *Tilting headpiece*
- *Adjustable ankle slide*
- *Variable height 21" to 29"*
- *Seamless body cushions 24" wide*
- *Hand rest and paper attachments*

Optional Features

- *Variable Height 18" to 28"*
- *Cervical manual drop-forward*
- *Cervical drop-toggle*
- *Thoracic Drop*
- *Pelvic Drop*
- *Elevating Headpiece (both cervical drops, lateral bending, lateral translation, foot activated elevation)*
- *Abdominal Breakaway*
- *Widths 24" or 26"*



TRU-DIGITAL

Imaging Solutions

1-844-944-6243

MIKE@TRU-DIGITAL-IMAGING.COM