



# Legacy 500

## Intersegmental Traction Table



### Table Specifications

Electrical: 120V, 60Hz,  
10amp

Weight Capacity 400lbs

Warranty 1 year limited

Lifetime warranty on  
rollers

Shipping Weight 349lbs

*With 2" more travel and almost 2" wider rollers than any one in the industry the Legacy-500 is specially designed to release muscle tension and deliver a soothing massage. Under the thick, soft, padded top, three large rotating rollers travel the back to release tension and relax tired muscles.*

### Standard Features

- 24" wide x 74" long
- 30-minute timer
- Reversible massage direction
- Adjustable massage pressure
- Vibration with On-Off feature
- Automatic massage pressure lowering
- Treat specific area or full spine
- Cervical roll and knee bolster

### Optional Features

- Reverse side controls
- Therapy top
- 27" over sized top
- Patient Straps and Rails
- Spring Loaded Roller Pressure
- Push Button Timer, 1-10 Minute Range
- Patient Hand Control - Stops Travel and control Roller Pressure



TRU-DIGITAL

Imaging Solutions

1-844-944-6243

MIKE@TRU-DIGITAL-IMAGING.COM