



SI- 1000



Table Specifications

Electrical: 120V, 60Hz,
10amp

Weight Capacity 400lbs

Warranty 1 year limited

Lifetime warranty on
rollers

Shipping Weight 349lbs

The SI-1000 has the same high quality as the Legacy-500 but is designed for the clinic where time is a factor. With an internal timer you never worry time will vary, a simple push button starts table to desired run time every time.

Standard Features

- 24" wide x 74" long
- 10-minute internal timer (set at desired run time)
- Push button start
- Counter-clock massage direction
- Adjustable massage pressure
- Automatic massage pressure lowering
- Treat full spine
- Cervical roll and knee bolster

Optional Features

- Reverse side controls
- Therapy top
- 27" over sized top
- Patient straps & rails
- Spring loaded roller pressure



TRU-DIGITAL

Imaging Solutions

1-844-944-6243

MIKE@TRU-DIGITAL-IMAGING.COM